


# μFOCAL VUV OPTICS

FOR UVS 300 AND TMM 304

## KEY FEATURES

- Retractable focusing capillary (< 250 μm/50 μm FWHM)
- Ultimate Photon Flux Densities
- UVS 300 μFOCAL high brilliance source (< 2 meV for He I)
- TMM 304 μFOCAL UV monochromator (< 1 meV for He I)



SPECS™

A member of SPECSGROUP

# μFOCAL VUV OPTICS

FOR UVS 300 AND TMM 304

## Focused VUV optics with 250 μm and 50 μm spot size for the proven high brilliance laboratory source UVS 300 and the UV monochromator TMM 304

Modern electron analyzers for Angle Resolved Photoelectron Spectroscopy (ARPES) experiments allow for highest angular resolutions. In complete experimental setup, the spot sizes of the laboratory UV sources used are meanwhile the limiting factor for further improvement. Additionally, the samples studied are often inhomogeneous at the surfaces with domains in the range of several hundreds of micrometers. An excitation of only a single domain also requires small VUV spot sizes. Furthermore, Momentum Microscopes benefit from small circular spots with adjustable sizes in the range between 50 μm and 500 μm, that can be adapted to the respective field of views.

The new series of μFOCAL 50 and μFOCAL 250 offers, as an addition to the proven range of focused VUV sources and monochromators, uncomparable and adjustable spot sizes. The use of special retractable capillaries makes them adaptable and easy to use to mostly all experimental

ARPES setups, allowing for circular spot sizes of well below 250 μm and even 50 μm at ultimate photon densities. By this, results taken in the laboratory are finally comparable in angular resolution and measuring times to experiments performed at modern beamlines/monochromators at synchrotron sources.




Fig 1. Spot Size of UVS 300/TMM 304 μFOCAL 250 measured in two perpendicular directions on calibrated photodiode for monochromated He I VUV light

### Technical Data

#### μFOCAL 250

Spot Diameter	Circular spot with 250 μm FWHM
Photon Flux	$> 3 \cdot 10^{11}$ photons/s
Photon Flux Density (monochromatized)	$> 1 \cdot 10^{14}$ photons/(s*mm) for He I, $> 6 \cdot 10^{12}$ photons/(s*mm) for He II
Photon line width	$< 2$ meV (He I radiation)
Working distance	10 mm
Operating pressure	$< 5 \cdot 10^{-8}$ mbar, $< 2 \cdot 10^{-10}$ mbar with monochromator
Monochromator	TMM 304 with two cassettes (line with $< 1$ meV for He I)

#### μFOCAL 50

Spot Diameter	down to 50 μm with up to 10 times higher flux density
---------------	---