


COMPONENTS FOR SURFACE ANALYSIS

Mini ECR Plasma Source **MPS-ECR**

- Microwave coupled plasma
- Based on 2.75" (NW40CF) flange
- For Oxygen, Nitrogen, Hydrogen
- Filamentless design
- No microwave tuning required
- User configurable coaxial design



The Mini ECR Plasma Source MPS-ECR

The SPECS Plasma Cracker Source MPS-ECR is a truly UHV compatible source for many PVD applications where a NW40CF base flange is necessary. The design is based on the well known SPECS PCS-ECR technology. A high density plasma is created by coupling a radially symmetric 2.45 GHz micro-wave field to ions on the 86 mT surface of a multi-polar magnetic array. A unique combination of features and options make this a highly versatile plasma source.

Main Features:

- **Mounting flange: 2.75" (NW40CF)**
- **UHV compatible**
- **Bakeable: 200° C**
- **Power: 60 W max. at 2.45 GHz**
- **Magnet type: permanent**
- **Integral water cooling (Feedthrough and jacket)**
- **In vacuum length: 130 mm**
- **Beam diameter: ~ 13 mm at source**
- **Gas flow rate: ~ 1 sccm typical**
- **Working distance: 50 - 300 mm**
- **Discharge tube: Materials according to gas**

Options:

- **Integrated shutter**
- **Various aperture types**
- **Differential pumping: Extends the range of working pressures into the ~ 10⁻⁷ and 10⁻⁸ mbar range**
- **Ion trap: Deflects residual ion current out of the beam.**

Operating Modes:

The source can be operated in three distinct modes, according to the extraction optics fitted and covering a wide range of ion energies and particle types.

Atom source – Thermal energy neutrals

This mode is intended for low energy and low damage surface treatment and sample growth. The particles released are largely thermalised (< 1 eV) and are therefore suitable for sensitive applications.

Downstream plasma / ECR mode – Low energy ions and neutrals

This mode reproduces classic ECR source characteristics. The optics here allow ions and higher energy plasma particles (~ 25 eV) to flood out into the chamber.

Hybrid – Atom source/ECR characteristics with controllable ions

This mode combines atom source and ion source behaviour to produce a source, which behaves like the atom source above until potentials are applied to the extraction grids, when ions with controllable energy (50 - 1500 eV) are then added to the beam.

Applications:

Nitrogen: Nitriding (GaN, AlN, InN and SiN), doping (ZnSe) and alloying (GalNAs, GaAlAsN)
Oxygen: HTc superconductors, optical coatings, dielectrics, reactive sputtering, laser ablation and ceramic growth (Al ₂ O ₃) oxygen cleaning and oxidation kinetics, post growth oxidation / low temperature SiO ₂
Hydrogen: Cleaning, growth enhancement / surfactant
Chlorine: In-situ etching
Methane (carbon): SiC film growth

SPECS GmbH - Surface Analysis
and Computer Technology
Voltastr. 5
13355 Berlin
Germany



Phone.: +49 (0)30 467824 - 0
Fax: +49 (0)30 4642083
E-mail: support@specs.de
http://www.specs.de

Your Representative: