


COMPONENTS FOR SURFACE ANALYSIS

Microwave Plasma Source **PCS-ECR**

- Microwave coupled plasma
- For Oxygen, Nitrogen, Hydrogen
- Filamentless design
- No microwave tuning required
- User configurable coaxial design
- Integral water cooling jacket



The Microwave Plasma Source PCS-ECR

The SPECS Plasma Cracker Source PCS-ECR is a truly UHV compatible source for most demanding MBE applications. Through the action of the well-known electron cyclotron resonance (ECR) phenomenon, a high density plasma is created by coupling a radially symmetric 2.45 GHz microwave field to ions on the 86 mT surface of a multi-polar magnetic array. The absence of a hot filament permits operation with most gases including reactive gases such as oxygen, chlorine, nitrogen and hydrogen. A unique combination of features and options make this a highly versatile plasma source.

Main Features:

- **Mounting flange: 4.5" (NW63CF)**
- **UHV compatible**
- **Bakeable: 200° C**
- **Power: 250 W max. at 2.45 GHz**
- **Magnet type: permanent**
- **Integral water cooling (feedthrough and jacket)**
- **In vacuum length: 300 mm**
- **Beam diameter: ~ 25 mm at source**
- **Gas flow rate: < 0.1 sccm to 100 sccm (mode dependant)**
- **Working distance: 50 - 300 mm**

Options:

- **Integrated shutter**
- **Non standard lengths**
- **Various aperture types**
- **Differential pumping: Extends the range of working pressures into the ~ 10⁻⁷ and 10⁻⁸ mbar range**
- **Faraday cup: Provides the possibility to monitor the beam current**
- **Ion trap: Deflects residual ion current out of the beam**

Operating Modes:

The source can be operated in three distinct modes, according to the extraction optics fitted and covering a wide range of ion energies and particle types.

Atom source – Thermal energy neutrals

This mode is intended for low energy and low damage surface treatment and sample growth. The particles released are largely thermalised (< 1 eV) and are therefore suitable for sensitive applications.

Downstream plasma/ECR mode – Low energy ions and neutrals

This mode reproduces classic ECR source characteristics. The optics here allow ions and higher energy plasma particles (~ 25 eV) to flood out into the chamber.

Hybrid – Atom source ECR characteristics with controllable ions

This mode combines atom source and ion source behaviour to produce a source, which behaves like the atom source above until potentials are applied to the extraction grids, when ions with controllable energy (50 - 1500 eV) are then added to the beam.

Applications:

Nitrogen: Nitriding (GaN, AlN, InN and SiN), doping (ZnSe) and alloying (GalnNAs, GaAlAsN)

Oxygen: HTc superconductors, optical coatings, dielectrics, reactive sputtering, laser ablation and ceramic growth (Al₂O₃) Oxygen cleaning and oxidation kinetics, post growth oxidation / low temperature SiO₂

Hydrogen: Cleaning, growth enhancement / surfactant

Chlorine: In-situ etching

Methane (carbon): SiC film growth

SPECS GmbH - Surface Analysis
and Computer Technology
Voltastr. 5
13355 Berlin
Germany



Phone.: +49 (0)30 467824 - 0
Fax: +49 (0)30 4642083
E-mail: support@specs.de
http://www.specs.de

Your Representative: